

2.0W LasePen



LasePen integrates the most state-of-art technology, delivers 2 watts of continuous laser power without any cords, outlets and footswitch. The laser communicates with the controller by a wireless modules. The control software is compatible with windows 7&XP operation system. You can easily install the software to your PC or an optional 7" Pad are available! 12 commonly used procedures have been programmed.

You can control the laser with simple one-touch presets, or set your own preferences, the choice is yours!

Applications

- Abscess
- Aphthous ulcer
- Biostimulation
- Bleaching
- Canal sterilization
- Contouring
- Curettage
- Frenectomy
- Gingivectomy
- Hemostasis
- Implant exposure
- Leukoplakia
- Sulcular Debridement
- Troughing

Technical Specifications

Lase Pen

Type of Laser:	Diode AlGaAs/InGaAsP
Wavelength:	810nm
Power Output:	0.1-2W
Aiming Beam:	650nm/1mW controllable
Operation Mode:	CW or Pulsed
Pulse Width:	0.5ms-30s
Dimension(LxD):	195mm xdia.18 mm
Weight:	<200g, with battery
Rechargeable Battery:	3.2 volt Lithium, 600mAh
Charging Station:	
Dimensions:	125 x70 x36 mm
Weight:	< 150g
Charging Ports:	4




Software available
or Optional 7" Pad



Features

- Compact, Pen style
- Finger Switch
- Battery-Operated
- Comprehensive Safety Precaution
- Full Accessories
- Ergonomic design

Big Role Small Size



Charging Station