

Model:PS800-A

SMOKE EVACUATOR

www.jhlaser.com

*Operate in
a forest.*



SHANGHAI JUE HUA LASER TECHNOLOGY DEVELOPMENT CO.,LTD

SMOKE EVACUATOR PS800-A



The PS800-A is a small, economical system designed to meet the needs of most medical specialties. It is light enough to easily move from room to room or be placed on a table.

Features

- Provides filtration is about 99.999%
- Extremely lightweight & portable
- Easy to operate
- Economical
- Simple filter changes
- Quiet

Technical specification

Product	Smoke evacuator
Model	PS800-A
Outlook	portable
Filter	4 stage ULPA filtering, VLSI level
Particle sizes	0.01
Efficient filtering	99.99995%
Motor power	Motor power: 1100W
Static evacuation mm. H ₂ O	1278
Footswitch control	Available
Display	LED display of filter condition, flow speed, footswitch control
Voltage	220V 50/60HZ
Size (cm ³)	18x40x30
Weight (Kg)	20

Consumable & accessory



Hose stand



Filter



Hose



2F, No.33, Lane2, Suide Rd.,
Shanghai 200331 P.R.China
Tel: +86 21 55214880
Fax: +86 21 55214080
<http://www.jhlaser.com>
E-mail: jhlaser@yahoo.com.cn